BETTER AG

Top quality at manufacturer prices

Bioteke COVID-19 Rapid test for professional use

With an integrated buffer solution

ID No.: 2067

Listed for EU-wide recognition in the "EU common list"
of the European Commission - Directorate General for Health and Food Safety
Common list of COVID-19 rapid antigen tests
[Click here](#) to check the validity of the CE certificate

Sensitivity	96,49%
Specificity	99,28%
Result in	15-20 Minutes
Packaging	25 pieces per box

It only takes 15 seconds for the antigen
to be released into the swab

Visible result in 15-20 Minutes

03/2022

09/2023

SARS-CoV-2 Antigen Test Kit (colloidal gold method)

User Instruction Manual

PRODUCT NAME

SARS-CoV-2 Antigen Test Kit (colloidal gold method)

PACKAGE SPECIFICATION

1 Test/Kit; 3 Tests/Kit; 5 Tests/Kit; 20 Tests/Kit; 50 Tests/Kit

INTENDED USE

This kit is only used for the *in vitro* qualitative detection of SARS-CoV-2 antigen from human nasopharyngeal swabs, oropharyngeal swabs and anterior nasal swabs specimens. SARS-CoV-2 Antigen Test Kit (colloidal gold method) is an immunochromatographic double-antibody sandwich assay intended for the qualitative detection of nucleocapsid protein antigen of SARS-CoV-2 from individuals who are suspected of COVID-19 by their healthcare provider within the first 7 days of symptom onset to suspect COVID-19 infection. This kit is suitable for the auxiliary diagnosis of COVID-19, the results are for clinical reference only and cannot be used as the sole basis for diagnosis and exclusion decision. The clinical diagnosis and treatment of patients should be considered in combination with their symptoms/signs, medical history, other laboratory tests and treatment responses. Positive test result needs to be further confirmed, negative result does not preclude SARS-CoV-2 infection. This kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of *in vitro* diagnostic procedures.

TEST PRINCIPLE

The kit is immunochromatographic and uses double-antibody sandwich method to detect SARS-CoV-2 antigen. During detection, the treated specimens are loaded into the sample wells of the test card. When the concentration of SARS-CoV-2 antigen in specimen is higher than the minimum detection limit, the viral antigen will form complexes with labeled antibodies first. Under chromatography, the complexes move forward along the nitrocellulose membrane till captured by pre-coated monoclonal antibody of SARS-CoV-2 in detection zone on nitrocellulose film (T) to form a pink/purple reaction line on the detection zone, at this point the result is positive; conversely, if there is no viral antigen or the concentration of antigen in specimen is below the minimum detection limit, no pink/purple reaction line appears in the detection zone, at this point the result is negative. Regardless of whether the sample contains viral antigens or not, a pink/purple reaction line will appear in the quality control zone (C), the pink/purple reaction line that appears in the quality control zone (C) is the criterion for determining if the chromatography process is normal.

MATERIALS PROVIDED

The test kit consists of test card, sample diluent, sample extraction tube, Waste bag, Oropharyngeal swab.

Components	Main Ingredients	Loading quantity (Specification)				
		1 Test/Kit	3 Tests/Kit	5 Tests/Kit	20 Tests/Kit	50 Tests/Kit
Test card	Test strip containing SARS-CoV-2 monoclonal antibody, Anti-mouse IgG polyclonal antibody	1 pc	3 pcs	5 pcs	20 pcs	50 pcs
Sample diluent		0.5mL	0.5 mL *3	0.5 mL *5	0.5mL*20	0.5mL*50
Sample extraction tube		1 pc	3 pcs	5 pcs	20 pcs	50 pcs
Waste bag		1 pc	3 pcs	5 pcs	20 pcs	50 pcs
Oropharyngeal swab		1 pc	3 pcs	5pcs	20 pcs	50 pcs

Note:

- Test cards are sealed together with desiccant in aluminum foil pouch.
- Do not mix use different batches of test cards and sample diluent.

STORAGE CONDITIONS AND SHELF LIFE

The test card and sample diluent should be stored at 2°C~30°C, valid for 18 months. Test cards should be used as soon as possible within 1 hour after opening the foil pouch. The bottle of sample diluent should be capped immediately after use and stored at 2°C~30°C, please use it within the validity period.

Date of manufacture and expiration: See package label for details.

SPECIMEN REQUIREMENTS

	1. Nasopharyngeal swab collection: Ask the patient to keep the head still to remove surface secretions from the anterior nasal foramen; Gently and slowly insert the swab through the nasal cavity to the nasopharynx, it reaches the posterior nasopharynx when resistance is encountered. Remain for several seconds to aspirate secretions; remove the swab with gentle rotation.
	2. Oropharyngeal swab collection: Insert the swab in the mouth completely into the pharynx, centering on the red swelling of the pharynx wall and upper anterior tonsils. Wipe both sides of pharyngeal tonsils and pharynx posterior wall with moderate force, avoid touching the tongue, and remove the swab.
	3. Anterior nasal swab collection: Ask the patient to keep the head still to remove surface secretions from the anterior nasal foramen. Insert the swab into one nostril of the patient. The swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nostril. 2. Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected. Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities. Withdraw the swab from the nasal cavity.

It is recommended that specimens be treated with the sample diluent provided with the kit as soon as possible after collection. If immediate processing is not possible, the specimen can be stored in a dry, sterilized and tightly sealed plastic tube at 2°C~8°C for up to 8 hours, and -70 °C for long-term storage.

TEST PROCEDURE

Please read the instruction manual seriously before testing. If the reagents in kit are stored in refrigerator, please take them out and leave them at room temperature before testing. The test should be done at room temperature.

I. Specimen extraction (see Figure 1).

- Twist the cap off
- Insert the sampled swab into the sample diluent in the sample extraction tube, squeeze the villi part of the swab in the tube through the outer wall of the tube by finger 5 times to dissolve the potential viral antigen into solution from swab as much as possible, then remove and discard the swab.
- Cover the dropper's lid on the sample extraction tube after step 2.

II. Testing procedures (see Figure 1).

- Remove the test card by opening the aluminum foil pouch from tear notch and lay it flat.
- Add 4 drops (approximately 80μL) of the treated specimen into the sample well of the test card.
- Please read the chromogenic results in the detection zone between 15~20 minutes to ensure proper test performance. Do not read results after 30 minutes. Results read after 30 minutes are invalid.


Figure 1

INTERPRETATION OF TEST RESULTS


Figure 2

- This kit contains quality control process, when the pink/purple reaction line appears in the C zone, it indicates the correct and effective operation. C line is the prerequisites to see if the test is valid. If the C line does not show color, regardless of whether the T line shows color or not, the test is invalid, and it is recommended to retest.
- The detection results of the kit are interpreted according to the following table.

C Line	T Line	Result
Pink/purple	Pink/purple	Positive
Pink/purple	Colorless	Negative
Colorless	Whether color is visible or not	Invalid test, retest

LIMITATIONS OF THE TEST

- The test results of this kit are only for the reference of clinicians and should not be used as the sole basis for clinical diagnosis and treatment. Clinical management of patients should be considered in the context of their symptoms/signs, medical history, other laboratory tests and response to treatment.
 - Sample collection and sample processing have a greater impact on the detection of pathogens, and a negative test result does not exclude the possibility of a viral infection.
 - Due to methodological limitations of antigen-based test, the analytical sensitivity of immunochromatographic method is generally lower than that of nucleic acid-based test. Therefore, the test provider should pay more attention to the negative results and make a comprehensive judgment based on other test results. It is suggested that the negative results in suspected patients should be checked by nucleic acid test or virus culture identification.
 - When the result of test kit is positive, it is recommended to combine the results of other methods (such as PCR and CT imaging) for further confirmation, and consult with local public health prevention institutions for treatment.
 - Analysis of the likelihood of false-negative results.
- (i) Improper sample collection, transport and processing, and low viral titers in the sample may lead to false negative results.
- (ii) The optimal sample type and the optimal sampling time after infection (peak viral titer) have not been validated, therefore, multiple sampling at multiple sites in the same patient may avoid false negatives.

PERFORMANCE CHARACTERISTICS

- The width of the membrane strip of this kit is not less than 2.5 mm, and the liquid migration speed is not less than 10 mm/min.
- Negative/positive reference coincidence rate

All the positive references are positive for the corresponding pathogens, which is consistent with the known results of the reference; all the negative references are negative for the corresponding pathogen.

3. Repeatability

Repeated testing was conducted for national or enterprise repeatable reference products for 10 times. The test results were consistent with the known results of the reference products and were uniform in color.

4. Limit of Detection (LoD)

The Limit of Detection (LoD) of SARS-CoV-2 Antigen Test Kit (colloidal gold method) is 1.75x10² TCID₅₀/mL.

5. Analytical specificity

- Sensitivity and Specificity

The SARS-CoV-2 Antigen Test Kit (colloidal gold method) was compared with the PCR Comparator Method on nasopharyngeal swab specimens.

Bioteke reagent	PCR reagent		Total
	Positive	Negative	
Positive	110	2	112
Negative	4	276	280
Total	114	278	392

Clinical sensitivity = 96.49% (95%CI: 91.26%~99.04%)

Clinical specificity = 99.28% (95%CI: 97.43%~99.91%)

Accuracy = 98.47% (95%CI: 96.70%~99.44%)

Kappa value = 0.9627

The SARS-CoV-2 Antigen Test Kit (colloidal gold method) Performance against PCR Comparator Method on oropharyngeal swab specimens.

BioTeke reagent	PCR reagent		Total
	Positive	Negative	
Positive	115	2	117
Negative	6	125	131
Total	121	127	248

Clinical sensitivity = A/(A+C)×100%= 95.04% (95%CI: 89.52%~98.16%)

Clinical specificity = D/(B+D)×100%= 98.43% (95%CI: 94.43%~99.81%)

Accuracy = (A+D)/(A+B+C+D)×100%= 96.77% (95%CI: 93.74%~98.60%)

Kappa value = 0.9354

The SARS-CoV-2 Antigen Test Kit (colloidal gold method) was compared with the PCR Comparator Method on anterior nasal swab specimens.

Bioteke reagent	PCR reagent		Total
	Positive	Negative	
Positive	131	1	132
Negative	12	163	175
Total	143	164	307

Clinical sensitivity = 91.61% (95% CI: 85.80%~95.59%)

Clinical specificity = 99.39% (95% CI: 96.65%~99.98%)

Accuracy = 95.77% (95% CI: 92.87%~97.73%); Kappa value = 0.9145

2) There is no cross-reactivity with the following pathogens: Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, MERS-coronavirus, SARS-coronavirus, Human coronavirus HKU1, Respiratory adenovirus type 1, Respiratory adenovirus type 2, Respiratory adenovirus type 3, Respiratory adenovirus type 4, Respiratory adenovirus type 5, Respiratory adenovirus type 7, Respiratory adenovirus type 55, Human Metapneumovirus (hMPV), Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Parainfluenza virus 4, Influenza A virus 2009H1N1, Influenza A virus seasonal H1N1, Influenza A virus H3N2, Influenza B virus Yamagata, Influenza B virus Victoria, Enterovirus 71, Respiratory syncytial virus Rhinovirus, Haemophilus influenza, Streptococcus pneumonia, Streptococcus pyogenes, Candida albicans, Pooled human nasal wash – representative of normal respiratory microbial flora, Bordetella pertussis, Mycoplasma pneumonia, Chlamydia pneumonia, Legionella pneumophila, Staphylococcus aureus, Staphylococcus epidermidis, Mycobacterium tuberculosis, Pneumocystis jirovecii (PJP), Measles virus, Human cytomegalovirus, Rotavirus, Norovirus, Mumps virus and Varicella zoster virus.

Pathogens likely to cross react with the test kit	Strain	Source/ Sample type
Human coronavirus 229E	ATCC VR-740, 229E	Virus cultures
Human coronavirus OC43	ATCC VR-1558, OC43	Virus cultures
Human coronavirus NL63	BELRESOURCES NR-470	Virus cultures
MERS-coronavirus	Synthesis	Pseudovirus
SARS-coronavirus	Synthesis	Pseudovirus
Human coronavirus HKU1	GZ/1804-138	Virus cultures
Respiratory adenovirus type 1	ADV1/GZ/Hecin1 608-21	Virus cultures
Respiratory adenovirus type 2	ADV1/GZ/Hecin1 609-42/2016	Virus cultures
Respiratory adenovirus type 3	ADV11/GZ/Hecin 1608-21	Virus cultures
Respiratory adenovirus type 4	ATCC VR-1572	Virus cultures
Respiratory adenovirus type 5	ADV3/GZ/1609-2	Virus cultures
Respiratory adenovirus type 7	ATCC VR-7	Virus cultures
Respiratory adenovirus type 55	ADV55/GZ/1612-129	Virus cultures
Human Metapneumovirus (hMPV)	GZ/1803-107	Virus cultures
Parainfluenza virus 1	PIV1/Guangzhou/0701/2011	Virus cultures
Parainfluenza virus 2	ATCC VR-92, Greer	Virus cultures
Parainfluenza virus 3	PIV3/Guangzhou/0902/2012	Virus cultures
Parainfluenza virus 4	ATCC VR-1378, M-25	Virus cultures
Influenza A virus 2009H1N1	A/GZ/GIRD02/2009(2009H1N1)	Virus cultures
Influenza A virus seasonal H1N1	A/PR/8/34(H1N1)	Virus cultures
Influenza A virus H3N2	A/Aichi/2/68(H3N2)	Virus cultures
Influenza B virus Yamagata	B/Guangzhou/GIRD06/09(Yamagata)	Virus cultures
Influenza B virus Victoria	B/Guangzhou/GIRD08/09(Victoria)	Virus cultures
Enterovirus 71	EV71/Guangzhou/0401/2011	Virus cultures
Respiratory syncytial virus	RSV A/GZ/Hecin1705-74	Virus cultures
Rhinovirus	A30/GZ/1710-89	Virus cultures
Haemophilus influenzae	GIM 1.961	Virus cultures
Streptococcus pneumoniae	GIM 1.550	Virus cultures
Streptococcus pyogenes	GDM1.546	Virus cultures
Candida albicans	FSCC 129002	Virus cultures
Pooled human nasal wash – representative of normal respiratory microbial flora	Clinical specimen	Clinical specimen
Bordetella pertussis	GDM 1.952	Virus cultures
Mycoplasma pneumoniae	ATCC 15531	Virus cultures
Chlamydia pneumoniae	ATCC VR-2282, TW-183	Virus cultures
Legionella pneumophila	GIM 1.840	Virus cultures
Staphylococcus aureus	FSCC 223005	Virus cultures
Staphylococcus epidermidis	GIM 1.143	Virus cultures

Mycobacterium tuberculosis	/	/
Pneumocystis jirovecii (PJP)	/	/
Measles virus	RA27/3	Virus cultures
Human cytomegalovirus	RC256	Virus cultures
Rotavirus	ATCC VR2018	Virus cultures
Norovirus	ATCC VR-3234SD	Virus cultures
Mumps virus	Jones	Virus cultures
Varicella zoster virus	CCTCC VR-1367	Virus cultures

3) Interfering substance: Human blood and mucins will not interfere with the results of the kit. The following common drugs will not interfere with the results of the kit: Oxymetazoline, Dexamethasone, Flunisolide, Sulphur, Kim Anh, Benzocaine, Zanamivir, Mupirocin, Tobramycin, Kalii Dehydrographolidi Succinas, Aspirin (enteric-coated tablets), Ibuprofen (granules), Acetaminophen (slow-release tablets) and Nimesulide Tablets.

Interfering substance		Test concentration
Mucins		1%
Human blood		5%
Nasal spray	Oxymetazoline	1.125mg/mL
Nasal corticosteroids	Dexamethasone	0.009mg/mL
	Flunisolide	0.75mg/mL
Zicam Cold Remedy Nasal Gel	Sulphur	335mg/mL
Allergic symptom relief drug	Kim Anh	4.5mg/mL
Oral anesthetic	Benzocaine	1.875mg/mL
Antiviral drug	Zanamivir,	75mg/mL
Antibiotics, nasal ointments	Mupirocin	33.5mg/mL
Systemic antibiotics	Tobramycin	0.3mg/mL
Immune system medication	Kalii Dehydrographolidi Succinas	3.8mg/mL
Antipyretic	Aspirin (enteric-coated tablets)	0.04g/L
	Ibuprofen (granules)	0.2g/L
	Acetaminophen (slow-release tablets)	450mg/L
	Nimesulide Tablets	0.05g/L

4) Hook effect: This kit doesn't have hook effect.

PRECAUTIONS

- 1. This is a single-use *in vitro* diagnostic reagent, do not reuse, and do not use expired products.
- 2. All test specimens must be considered potentially infectious, and during collection, processing, storage, mixing of specimens and testing should be taken appropriate protective measures. Therefore, wear protective measures such as wearing gloves and masks should be done, and dispose of all wastes as potentially biohazardous items. (Used cotton swabs, test cards, extraction tubes, etc., should be decontaminated before disposal and tested for autoclaving.)
- 3. Use the swab and sample diluent provided with this reagent for sampling, and do not mix use different batches of test cards and sample diluent.
- 4. Use fresh specimens for testing, do not use repeated freeze-thaw samples.
- 5. Operate at room temperature. Test cards kept at low temperature should be restored to room temperature before opening to avoid moisture absorption.
- 6. Do not use reagent kits with obvious damage or test cards with damaged or expired packaging.
- 7. The aluminum foil pouch contains desiccant and must not be taken orally.
- 8. Improper sample collection or processing may result in false-negative results.
- 9. Ensure proper sample loading volume, results of too much or too little sample loading volume may not be credible.
- 10. If the initial screen is a positive sample, contact your local public health agency.
- 11. As with the diagnostic reagents used, the final diagnosis should be made by a physician after combining the various test parameters and clinical symptoms.
- 12. If you have any questions or suggestions on the use of this kit, please contact the manufacturer.
- 13. For unknown reasons, long-term use of some drugs may lead to false positive results of the test, which are not covered by the interfering substances.

SYMBOLS

Date of manufacture


Do not use if package is damaged


CE mark of conformity

Keep away from sunlight

For *In Vitro* Diagnostic Use


Batch code

Keep dry

Consult instructions for use


Temperature limit

Do not re-use

Authorized representative in the European Community

Use-by date

Manufacturer

BIOTEKE CORPORATION (WUXI) CO., LTD
4th Floor, D5&2nd Floor, D3 & 1st and 2nd Floor, D16, No.1719,
Huishan Avenue, Wuxi, Jiangsu, CN 214174.

 SUNGO Europe B.V.
Olympisch Stadion 24,1076 DE Amsterdam ,Netherlands
Importer: Better AG
General-Guisan-Str. 8
6300 Zug, Switzerland
Tel: + 353 1 513 7511
Email: info@OdemShop.com
Shop: www.OdemShop.com

Bioteke Corporation(wuxi) Co.,Ltd
Address: 4th floor, D5, No.1719,Huishan Avenue,Wuxi,China

Manufacturer's Declaration

To whom it may concerns,

Product name: SARS-CoV-2 Antigen Test Kit (colloidal gold method)

Country of Origin: China

We, Bioteke Corporation(wuxi) Co., Ltd, headquartered in, 4th floor, D5, No.1719, Huishan Avenue, Wuxi, China, do hereby declare "Better AG" located in General-Guisan-Str. 8, 6300 Zug, Switzerland, is authorized to import, sell, distribute the "Bioteke" branded in Europe and Africa.

We hereby confirm the authenticity of the SARS-CoV-2 Antigen Test Kit (colloidal gold method) sold by this distributor

Bioteke Corporation(wuxi) Co., Ltd

June 28, 2022

**Effectiveness Statement of BioTeke Novel Coronavirus (COVID-19)
test kit for the SARS-CoV-2 Variants Detection**

In quick response to the newly found variant of novel coronavirus Omicron (B.1.1.529) in Botswana, we have analyzed our Novel Coronavirus (COVID-19) test kit and the result shows that:

Our Freeze-dried Novel Coronavirus (COVID-19) nucleic acid detection kit (PR2019, PR2019-D, and PR2020-D) and SARS-CoV-2 Antigen Test Kit (colloidal gold method) (TC1002) can detect the new variant B.1.1.529 accurately, and for variants, which include but not limited to Alpha (B.1.1.7), Beta (B.1.351), Gamma (P.1), Delta (B.1.617.2), Lambda (C.37), Kappa (B.1.617.1), Eta (B.1.525), Lota (B.1.526), Mu (B.1.621), Zeta (P.2), and Omicron (B.1.1.529) etc., the results did not show any off-target and missing detection. The accuracy and sensitivity can be guaranteed.

BioTeke will continue to track and quickly respond to the latest variant of novel coronavirus and ensure that there will be no off-target and missing detection of our test kit.

BioTeke Corporation(wuxi) Co.,Ltd

November 27, 2021

HZJZ

HRVATSKI ZAVOD

ZA JAVNO ZDRAVSTVO

Rockefellerova 7

HR-10000 Zagreb

T: +385 1 4863 222

F: +385 1 4863 366

www.hzjz.hr

VALIDATION STUDY REPORT

OBJECTIVE:

To establish the performance (sensitivity and specificity) of BioTeke SARS-CoV-2 Antigen Test Kit, BIOTEKE CORPORATION (WUXI) CO.,LTD, China (LOT No.: TC1002E2022L0404; MDD: 02/2024).

Samples included:

100 RT-PCR SARS-CoV-2 positive samples from persons within the first 7 days from the onset of COVID-19 symptoms during the Omicron wave (November 19-24, 2022)

50 RT-PCR SARS-CoV-2 negative symptomatic or asymptomatic persons.

Samples were taken at the drive-in test station of the Croatian Institute of Public Health from symptomatic and asymptomatic persons. The samples were collected by healthcare professionals. Nasopharyngeal swabs were taken and placed into the viral transport medium (2 ml). The swabs were transported to the laboratory and tested with the COVID-19 Antigen Rapid Test and with confirmatory RT-PCR test (Allplex™ 2019-nCoV Assay, Seegene Inc, Republic of Korea, Lot:RV9121L95, Exp.:1/2023). RT-PCR and COVID-19 Antigen Rapid Test were performed from the same nasopharyngeal swab.

Diagnostic acceptance criteria:

- Assay should have sensitivity over 80%, or over 90% for subjects with a Ct < 25
- Assays should have a specificity over 98%.

Ethical approval

The study was conducted according to the guidelines of the Declaration of Helsinki and approved by the Ethics Committee of the Croatian Institute of Public Health (Protocol No. 030-02/22-05/4, Approved on October 28, 2022).

Results:

	POSITIVE RT-PCR (all, Ct<30; Ct<25)	NEGATIVE RT-PCR	TOTAL
POSITIVE RAT	67/67/65	0	67
NEGATIVE RAT	33/17/6	50	83/67/56
TOTAL	100/84/71	50	150/134/121

All samples

	Statistic Value	95% CI
Sensitivity	67.00%	56.88% to 76.08%
Specificity	100.00%	92.89% to 100.00%
Negative Likelihood Ratio	0.33	0.25 to 0.44
Positive Predictive Value (*)	100.00%	
Negative Predictive Value (*)	60.24%	53.40% to 66.70%
Accuracy (*)	78.00%	70.51% to 84.35%

Samples (Ct<30)

	Statistic Value	95% CI
Sensitivity	79.76%	69.59% to 87.75%
Specificity	100.00%	92.89% to 100.00%
Negative Likelihood Ratio	0.20	0.13 to 0.31
Positive Predictive Value (*)	100.00%	
Negative Predictive Value (*)	74.63%	65.80% to 81.81%
Accuracy (*)	87.31%	80.47% to 92.43%

Samples (Ct<25)

	Statistic Value	95% CI
Sensitivity	91.55%	82.51% to 96.84%
Specificity	100.00%	92.89% to 100.00%
Positive Predictive Value (*)	100.00%	
Negative Predictive Value (*)	89.29%	79.49% to 94.71%
Accuracy	95.04%	89.52% to 98.16%

The CT values for the tested samples were ≤ 34 . Statistical analysis was done using online tool (https://www.medcalc.org/calc/diagnostic_test.php).

The assay has sensitivity over 90% for subjects with a Ct < 25 (91.55%) and specificity 100.00% which is in line with the requirements MDCG Guidance on performance evaluation of SARS-CoV-2 *in vitro* diagnostic medical devices.

Zagreb, January 05, 2023

Irena Tabain, PhD, MD
Head of Virology Department
Croatian Institute of Public Health

CE Certification – CIBG Registration Letter

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
Dr. M.J. van de Velde

CE Certification-EC Declaration of Conformity

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DECLARATION OF CONFORMITY

Regarding In Vitro Diagnostic Directive (98/79/EC)

Manufacturer: BioTeke Corporation (Wuxi) Co., Ltd

Address: 4th Floor, D5&2nd Floor, D3& 1st and 2nd Floor, D16, No.1719, Huishan Avenue,
Wuxi, JiangSu, CN 214174.

EC Representative: SUNGO Europe B.V.

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Product Name: SARS-CoV-2 Antigen Test Kit (colloidal gold method)

Specification: TC1002 (1 Test/Kit; 20 Tests/Kit; 50 Tests/Kit)

Classification: Others (IVDD)

Conformity Assessment

Procedure: Annex III of In Vitro Diagnostic Directive (98/79/EC)

We herewith declare that the above-mentioned products meet the requirements of In
Vitro Diagnostic Directive (98/79/EC) and the following harmonized standards.

EN ISO 14971:2012	EN ISO 18113-1:2011	EN ISO 18113-2:2011
EN 13612:2002+AC:2002	EN ISO 23640:2015	EN 13641:2002

Signature: 
Name/ Position: Zhitu Zhou / GM

Date: 9th November, 2020

Place: Wuxi, Jiangsu / China

On behalf of SUNGO Europe office, I confirmed we are
EU REP of the company who issue this document.

Authorized Signature (S)

MHRA Registration- UKCA Declaration of Conformity

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UKCA	
DECLARATION OF CONFORMITY	
Regarding UK Medical Devices Regulations 2002	
Manufacturer:	BioTeke Corporation (Wuxi) Co., Ltd
Address:	4 th Floor, D5&2 nd Floor, D3 & 1 st and 2 nd Floor, D16, No.1719, Huishan Avenue, Wuxi, Jiangsu, CN 214174
UK Responsible Person:	SUNGO Certification Company Limited
Address:	3rd floor, 70 Gracechurch Street, London. EC3V 0HR
Product Name:	SARS-CoV-2 Antigen Test Kit (colloidal gold method)
Model:	TC1002 (1 Test/Kit; 20 Tests/Kit; 50 Tests/Kit)
Category:	General IVDs
Conformity Assessment Procedure:	Annex III, Part IV of the UK MDR 2002
We herewith declare that the above-mentioned products meet the requirements of UK Medical Devices Regulations 2002 and the following standards.	
BS EN ISO 14971:2019	BS EN ISO 23640: 2015
BS EN ISO 15223-1:2016	BS EN ISO 18113-1: 2011
BS EN 13612: 2002/ AC:2002	BS EN ISO 18113-2: 2011
BS EN 13641: 2002	
Signature: 	Date: 2021/7/1
Name / Position: Zhitu Zhou / GM	Place: Wuxi, Jiangsu / China
	
The declaration of conformity is valid in connection with the release technical document BTK-IVD/UKCA-TCF001.	
All the supporting documentation is retained at the premises of the manufacturer. The Declaration of Conformity is exclusively under the sole responsibility of the manufacturer.	

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* POC = Point of Care

11.06.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with CT<25, 23 samples with CT between 25 and 30, and 9 samples with CT>30. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 µL of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

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EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY
Public health, country knowledge, crisis management
Health Security

EU health preparedness:

A common list of COVID-19 rapid antigen tests and a common standardised set of data to be included in COVID-19 test result certificates

Agreed by the Health Security Committee

This document was agreed by the HSC on 17 February 2021

Annex I

Common list of COVID-19 rapid antigen tests

A first update was agreed by the HSC on 10 May 2021; A second update was agreed by the HSC on 16 June 2021; A third update was agreed by the HSC on 7 July 2021; A fourth update was agreed by the HSC on 14 July 2021

IMPORTANT: A (interim) grace period of 8 weeks applies whenever updates are made to Annex I, the common list of COVID-19 rapid antigen tests

Annex II

Common standardised data set of to be included in COVID-19 test result certificates

An update to Annex II was agreed by the HSC on 19 March 2021

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Manufacturer	RAT commercial name	CE marking	Clinical performance <i>Data by manufacturer</i>	Clinical performance <i>Data used in MS</i>	FIND evaluation studies	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	In JRC database (Device ID #) ¹¹
				DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100% NL: Independent field study, mainly symptomatic individuals, sensitivity Ct≤30: 96.0%; specificity overall: 100%						
BIOSYNEX SA	BIOSYNEX COVID-19 Ag+ BSS	Yes	Clinical Sensitivity: 97.5 %	FR: Validation study data: 125 positive and 118 negative samples; sensitivity 96%, specificity: 99%		FR		FR		Yes (1494)
BIOTEKE CORPORATION (WUXI) CO., LTD	SARS-CoV-2 Antigen Test Kit (colloidal gold method)	Yes	96.49 % sensitivity 99.28 % specificity OP/NP swab	DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct25) + Manufacturer specificity: 99.28%		DE ^[2]		DE ^[2]		Yes (7067)
BTNX Inc	Rapid Response COVID-19 Antigen Rapid Test	Yes	90.2% sensitivity 100% specificity NP swab, NP swab, OP swab	DE: 94.55% sensitivity, 100% specificity		AT, DE ^[2] , ES, SI		DE ^[2]		Yes (1236)
CerTest Biotec	CerTest SARS-CoV-2 Card test	Yes	92.9% sensitivity 99.6% specificity NP swab	ES: Ct ≤ 25, sensitivity: 94,0%; sensitivity for samples within the first 5 days after symptom onset: 84,8%		ES, PT, SI		DE ^[2] , ES		Yes (1173)
Core Technology Co., Ltd	Coretests COVID-19 Ag Test	Yes	98.1% sensitivity 99.6% specificity NP swab	DE: 98.1% sensitivity, 99.6% specificity		AT, DE ^[2] , RO		DE ^[2]		Yes (1919)
CTK Biotech, Inc	OnSite COVID-19 Ag Rapid Test	Yes	Clinical Sensitivity: 92.3 % Clinical Specificity: 100 % Nasal, NP swab	ES: 219 samples; Nasal swab - Clinical sensitivity 86% (90%: Ct <30) Specificity: 100% (Method B) DK: 107 samples; Nasal swab - clinical sensitivity 86%; (from asymptomatic and mild symptomatic individuals), Clinical specificity: 100%	To start	DK		DK, ES		Yes (1581)
DDS DIAGNOSTIC	Test Rapid Covid-19 Antigen (tampon nazofaringian)	Yes	98.77% sensitivity 99.03% specificity Nasal swab	RO: Meets the minimum performance requirements.		RO		RO China	RO	Yes (1225)

Czech List of antigenic tests for which issued by the ministry exemption 捷克白名单

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SARS-CoV-2 Antigen Rapid Test Kit (Fluorescence Immunochromatography)	BIOHIT HealthCcare (Hefei) Co., Ltd.
SARS-CoV-2 Ag Rapid Test	BioMaxima SA
Biomerica COVID-19 Antigen Rapid Test (nasopharyngeal swab)	Biomerica Inc.
NowCheck COVID-19 Ag Test	BIONOTE
CORONAVIRUS AG RAPID TEST CASSETTE	BIO-RAD
BIOSYNEX COVID-19 Ag BSS	BIOSYNEX S.A.
BIOSYNEX COVID-19 Ag+ BSS	BIOSYNEX SA
SARS-CoV-2 Antigen Test Kit (colloidal gold method)	BIOTEKE CORPORATION (WUXI) CO., LTD
biotical SARS-CoV-2 Ag Card	Biotical Health S.L.U.BIOTICAL HEALTH S.L.U
AFLAS COVID-19 Ag	Boditech Med Inc
Rapid Response COVID-19 Antigen Rapid Test	BTNX Inc
CerTest SARS-CoV-2 Card test	CerTest Biotec
Coretests COVID-19 Ag Test	Core Technology Co., Ltd
OnSite COVID-19 Ag Rapid Test	CTK Biotech, Inc
Test Rapid Covid-19 Antigen (tampon nazofaringian)	DDS DIAGNOSTIC
DIAQUICK COVID -19 Ag Cassette	DIALAB GmbH
COVID-19 Antigen Detection Kit	DNA Diagnostic
Edinburgh Genetics ActivXpress+ COVID-19 Antigen Complete Testing Kit	Edinburgh Genetics Limited
EBS SARS-CoV-2 Ag Rapid Test	Eurobio Scientific
ESPLINE SARS-CoV-2	Fujirebio
GA CoV-2 Antigen Rapid Test	GA Generic Assays GmbH
Genbody COVID-19 Ag Test	GenBody Inc
SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	Genrui Biotech Inc
GenSure COVID-19 Antigen Rapid Test Kit	GenSure Biotech Inc
SARS-CoV-2 Antigen (Colloidal Gold)	Getein Biotech, Inc
One Step Test for SARSCoV-2 Antigen (Colloidal Gold)	Getein Biotech, Inc.
SARS-CoV-2 Antigen Kit (Colloidal Gold)	Goldsite Diagnostic Inc.
GENEDIA W COVID-19 Ag	Green Cross Medical Science Corp.
2019-nCoV Antigen Test Kit (colloidal gold method)	Guangdong Hecin Scientific, Inc.
COVID-2019-nCoV Ag Rapid TestDetection Kit (ImmunoChromatography)	Guangdong Longsee Biomedical Co., Ltd.
COVID-19 Ag Test Kit	Guangdong Wesail Biotech Co. Ltd

France COVID-19 Diagnostic test list - For nasopharyngeal swab

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<https://covid-19.sante.gouv.fr/tests>

Belgium famhp List of the recommended tests

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nl fr en Other information and services: www.belgium.be .be

famhp federal agency for medicines and health products

About the FAMHP Jobs Publications Press Contact Complaints Webportal

Your medicines and health products, our concern

Search

Home Human use Veterinary use Information for the public Information for professionals

Home ▶ Human use ▶ Health Products ▶ Medical devices and their accessories ▶ COVID-19 ▶ Tests

Tests

Information for actors (manufacturers, importers, recorded distributors)

- [Obligation to register antigenic rapid tests and antigen self-tests made available in Belgium by manufacturers, importers, recorded distributors](#)

Antibody and antigen tests for professional use

- [New validation procedure for serological tests and antigenic tests \(12/08/2021\)](#)
- [List of the recommended tests](#)

Self-tests

- [Validation procedure for antigen self-test with a CE-certificate](#)
- [List of recommended SARS-CoV-2 antigen self-tests with a CE-certificate](#)
- [List of professional rapid tests that may be sold as self-tests](#)

Nose swabs

- [Swabs: compliance verification](#)

Last updated on 17/12/2021

Notification of adverse reactions or incidents

PIL and SPC of a medicine

(PIL: patient information leaflet SPC: summary of product characteristics)

Search

News

21/12/2021
Limited availability of subcutaneous immunoglobulins

09/12/2021
Recall of ultrasound gel from Eco-Med Pharmaceuticals Inc.

99	Bioteke Corporation (Wuxi) Co	SARS-CoV-2 Antigen Test Kit (Colloidal gold method)	96.5	99.3	NP swab	1463
2067						

https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/covid_19/tests

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Holds Certificate No:

MDSAP 751219

Statement of Conformity: The company listed on this certificate has been audited to and found to conform with the following criteria: ISO 13485:2016 and Australia - Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure; Brasil - RDC ANVISA n. 16/2013, RDC ANVISA n. 23/2012, RDC ANVISA n. 67/2009; Canada - Medical Devices Regulations - Part 1 - SOR 98/282; Japan - MHLW Ministerial Ordinance 169, Article 4 to Article 68, PMD Act; USA - 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D

Design and Development, Production, Distribution and Service of Automated Nucleic Acid Extraction System, Nucleic Acid Extraction Kits, Disposable virus sampling Swab kits, Colloidal Gold test kits, Immunoassay Test Kits, Fluorescence PCR test kits and PCR in-vitro diagnostic instrument.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2021-11-28

Effective Date: 2021-11-28

Expiry Date: 2024-11-24

BSI Group America Inc. is an MDSAP authorized auditing organization

Page: 1 of 2

...making excellence a habit.™

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

ISO 13485 Certificate-MDSAP

ISO 13485 证书-MDSAP

Certificate No: **MDSAP 751219**

Location	Registered Activities
BioTeke Corporation(wuxi)Co., Ltd 4th Floor D5 No.1719 Huishan Avenue Wuxi City, Jiangsu P.R 214174 China Facility ID Number: F005548	Design and Development, Production, Distribution and Service of Automated Nucleic Acid Extraction System, Nucleic Acid Extraction Kits, Disposable virus sampling Swab kits, Colloidal Gold test kits, Immunoassay Test Kits, Fluorescence PCR test kits and PCR in-vitro diagnostic instrument.
BioTeke Corporation(wuxi)Co.,Ltd 2nd Floor, D3 No.1719,Huishan Avenue Wuxi City China Facility ID Number: F005548	Manufacture for Nucleic Acid Extraction kits.
BioTeke Corporation(wuxi)Co.,Ltd 1st and 2nd Floor D16,No.1719,Huishan Avenue WuXi Jiangsu 214174 China Facility ID Number: F005548	Manufacture for Automated Nucleic Acid Extraction System and Disposable virus sampling Swab kits.
BioTeke Corporation(wuxi)Co.,Ltd No. 330, Qiyang South Road, Jiangyin Qingyang Jiangsu 214401 China Facility ID Number: F005548	Manufacture for SARS-CoV-2 Antigen Test Kit (colloidal gold method), Nucleic Acid Extraction Kit ,Freeze-dried Novel Coronavirus (COVID-19) Nucleic Acid Detection Kit (Fluorescence PCR) ,Novel Coronavirus (SARS-CoV-2) nucleic acid detection kit Fluorescence PCR method, and Disposable virus sampling Swab kits.

Original Registration Date: 2021-11-28

Effective Date: 2021-11-28

Expiry Date: 2024-11-24

Page: 2 of 2

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](http://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

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CERTIFICATE

The Governing Board of
ARES International Certification Co., Ltd.
Hereby Grants To:

BIOTEKE CORPORATION (WUXI) CO., LTD.

Organization Credit Code: 913202065617502076

4th Floor, D5&2nd Floor, D3& 1st and 2nd Floor,D16, No.1719, Huishan Avenue,
Wuxi

Has been assessed and found to be in accordance with the requirements of standard
detailed below

GB/T 45001-2020/ISO 45001:2018

Scope

**Manufacture of disposable virus sampling swab and nucleic acid
extraction Management of Related Occupational Health and Safety
Aspects.**

Certificate No.: ARES/CN/12101044S

Certificate Issue Date: 2021-02-07

Registration Expiration Date: 2024-02-06

The time interval between each surveillance audit and the last on-site audit shall not exceed 12 months, and the
organization must obtain "surveillance audit approval notification" issued by ARES to ensure the validity of the certificate.

Authorized by :

Chiora Juwen

ARES International Certification Co., Ltd.

No.12-2, Ln. 187, Wenping Rd., Anping Dist., Tainan City 708, Taiwan

TEL / 06-295 9696 (Rep. Line) FAX / 06-295 9667 www.ares-registration.com

Check the validity of this certificate on the official website of Certification and Accreditation Administration
of the People's Republic of China (www.cnca.gov.cn) or www.ares-china.cn.

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ENVIRONMENTAL MANAGEMENT SYSTEM CERTIFICATE

Certificate No.: GXZT002-21E10015R0S

We hereby certify that the organization:

Bioteke (wuxi) Corporation Co., Ltd

Unified social credit code/Organization code: 913202065617502076

is in conformity with Environmental Management System Standard:

GB/T24001-2016/ISO14001:2015

The certificate is valid to the following product(s)/service:
the manufacture of disposable virus sampling swab kits and nucleic acid extraction
system and related management action

Registration Address: 4th floor, D5, 3rd floor, D3, 1st & 2nd floor, D16, No.1719,
Huishan Avenue, Wuxi

Audit Address: 4th floor, D5, 3rd floor, D3, 1st & 2nd floor, D16, No.1719,
Huishan Avenue, Wuxi

Date of Issue : 14-01-2021

Date of Expiry: 13-01-2024

The Effectiveness of the Certificate is Subject to QR Code in the Left.
Meanwhile, You Can Search the CNCA Website: www.cnca.gov.cn
or Website of Certification Body www.isogx.cn

Guo Xin Zheng Tong (Beijing) Inspection & Certification Co., Ltd.

Room 508, Building 42, Zone 2, Tiantongyuan, Changping District,
Beijing, China (102218)

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Further references and information

BETTER AG
General-Guisan-Str.8
6300 Zug, Switzerland

SCAN ME

IRL: +353 1 513 75 11
CH: + 41 71 58 80 248
Shop: www.OdemShop.com
E-Mail: info@OdemShop.com