BETTER AG


Top quality at manufacturer prices

COVID-19 quick test for professional use

With integrated buffer solution

Listed for EU-wide acceptance in the "EU-common list" of the European Commission - Directorate General for Health and Food Safety Common List of COVID-19 Antigen Rapid Tests

The result is visible after 15-20 minutes.

Sensitivity	96.77%
Specificity	100%
Result after	15 - 20 Minutes
Packing	25 pcs per box

Shenzhen Lvshiyuan Biotechnology Co.,Ltd

Green Spring® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Instructions for Use

REF GF102B1L Rev. 7.1

English

Rapid test for the qualitative detection of SARS-CoV-2 nucleocapsid antigen. For professional use.

INTENDED USE

The Green Spring® SARS-CoV-2 Antigen Rapid Test is intended for the rapid qualitative detection of SARS-CoV-2 nucleocapsid protein antigen in human saliva, nasal, nasopharyngeal or oropharyngeal swab samples. The results are used for the detection of SARS-CoV-2 antigen. The antigen is generally detectable in upper respiratory tract samples during the acute phase of infections. Positive results do not exclude bacterial infection or co-infection with other viruses. The pathogen detected may not be the sole cause of the disease.

Negative results should be treated as suspected cases and confirmed with a molecular assay. Negative results should be considered in the context of a patient's recent exposures, history and presence of clinical signs and symptoms consistent with COVID-19. The test should only be performed by trained medical personnel.

SUMMARY

The novel coronaviruses belong to a β -genus. COVID-19 is an acute respiratory infectious disease. Humans are generally susceptible to it. Currently, patients infected with the novel coronavirus are the main source of infection; asymptotically infected people may also be a source of infection. The main manifestations include fever, fatigue and a dry cough. A stuffy or runny nose, sore throat, muscle aches and diarrhoea occur in a few cases.

TEST PRINCIPLE

The Green Spring® SARS-CoV-2 Antigen Rapid Test is a qualitative, membrane-based immunoassay for the detection of SARS-CoV-2 nucleocapsid protein antigen. The test line area is coated with SARS-CoV-2 antibody. The sample reacts with the SARS-CoV-2 antibody in the test line area. If the sample contains SARS-CoV-2 antigen, a coloured line appears in the test line area as a relevant result. As a procedural control, a coloured line appears in the control line area, indicating that the correct volume of sample has been applied and membrane wetting has proceeded correctly.

STORAGE AND STABILITY

Store the tests in the sealed foil pouch at room temperature or refrigerated (2 - 30 °C). The test is stable until the expiry date printed on it. The test cassettes must be stored in the sealed foil pouch until use. Do not freeze. Do not use after the expiry date. Protect from sun, moisture and heat.

MATERIALS SUPPLIED

- Test cassettes (25 pieces)
- Sampling swabs: 25 pieces

- Extraction tubes with buffer: 25 pieces disposable reaction tubes with 0.5 ml extraction buffer each and 25 pieces nozzle cap
- Package leaflet: 1 instruction leaflet
- Workstation: 1 piece
- Desiccant: 1 package

PRECAUTIONS

1. Read the package leaflet carefully before performing the test. Failure to follow the instructions in the package leaflet may result in inaccurate test results.
2. For professional in vitro diagnostic use only. Do not use after the expiry date.
3. Do not eat, drink or smoke for 10 minutes before and during sample collection.
4. Do not use the test if the packaging or test components are damaged.
5. All samples must be considered potentially infectious. Observe established precautions against microbiological hazards throughout the collection, handling, storage and disposal of patient samples and used test components.
6. Wear protective clothing such as lab coats, disposable gloves and eye protection while testing samples.
7. Wash your hands thoroughly after performing the test.
8. Viral transport media (VTM) may affect the test result: Extracted samples for PCR testing cannot be used for testing.
9. All used test components should be disposed of according to local regulations.
10. Humidity and temperature may adversely affect the results.

PREPARATION

Use only the materials supplied with the respective set. Test the samples immediately.

Use the test kit only at room temperature (15 to 30 °C). The test kit is intended only for swab samples that are collected and tested directly (i.e. swabs that have NOT been placed in transport media). This kit is NOT intended for testing liquid samples such as wash or aspirate samples or swabs in transport media, as results may be affected by over-dilution.

1. Tear off the foil pouch, take out the test cassette and place it on a clean and flat surface.
2. Freshly collected samples should be processed within 1 hour.
3. Label the respective test cassette for each test or control.
4. Place the labelled extraction tubes in a rack in the designated area of the workspace.

COLLECTING THE SAMPLE

Correct sample collection is the most important step. Select one of the four methods and then proceed with the test procedure.

1) Saliva (lollipop)

Be aware that false results may occur if the saliva is not collected properly.

1. Place an extraction tube in the cardboard workstation.

2. Press the tip of your tongue against the lower root of your jaw. Cough deeply. Make the sound of "kuuuu" to gather the saliva.
3. Place the swab on the tongue for at least 10 seconds, rotating it 3 times or more to fully absorb the saliva.

2) Anterio-nasal swab (nose in front).

Make sure to collect enough nasal secretions with the swab. It is advisable to blow your nose first.

1. Place an extraction tube in the cardboard workstation.
2. Carefully insert the swab into the patient's nostril. The tip of the swab should be inserted up to 2.5 cm deep from the edge of the nostril.
3. Swab along the mucosa in the nostril to ensure that both mucus and cells are collected.
4. Take the swab out of the nostril while gently rotating it between your fingers.

3) Nasopharyngeal swab (nose-throat).

1. Place an extraction tube in the cardboard workstation.
2. Tilt the patient's head slightly backwards. Hold the swab like a pen and insert it through the nostril parallel to the palate.
3. While inserting, gently rub and roll the swab. As soon as you feel the throat resistance, stop and let the swab absorb secretion.
4. Slowly and carefully move the swab outwards while gently rotating it between your fingers.

4) Oropharyngeal swab (throat).

1. Place an extraction tube in the cardboard workstation.
2. Let the patient open the mouth wide and make "Ah" sounds, which will expose the pharyngeal tonsils on both sides.
3. Hold the swab tightly and wipe back and forth on the pharyngeal tonsils on both sides at least three times per side with moderate force. Do not touch the palate, tongue, teeth or gums.
4. Remove the swab while gently rotating it between your fingers.

For best results, the nasopharyngeal (nose-throat) method is recommended.

PERFORMING THE TEST

After collecting the sample, perform the test as follows:

1. Tear off the sealing of the extraction tube.
2. Insert the swab sample into the extraction tube and dip it up and down in the liquid. Rotate the swab several times during this process.
3. While removing the swab squeeze the sides of the tube to extract the remaining liquid out of the swab. Place the dropper tip firmly on the extraction tube and mix the liquid thoroughly.
4. Dispense 3 drops (approximately 100uL) into the sample well of the test cassette via the dropper tip.
5. Interpret the test results after 15 minutes. Do not interpret the results after 20 minutes.

INTERPRETING THE TEST RESULT

POSITIVE: Two lines appear. One coloured line appears in the control line area (C) and another coloured line appears in the test line area (T). A positive result in the test area indicates the detection of SARS-CoV-2 antigen in the sample. A positive result does not exclude infection with other pathogens.

NEGATIVE: A coloured line appears in the control area (C). No coloured line appears in the test line area (T). A negative result does not exclude viral infection with SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.

INVALID: The control line does not appear. Insufficient sample volume or incorrect handling are the most likely reasons causing the

control line not to appear. Check the procedure and repeat the test with a new test cassette. If the problem persists, stop using the test kit immediately and contact your dealer.

QUALITY CONTROL

The control area (C) acts as an internal procedure control. A coloured line appears when the procedure or sample volume has been applied correctly. Control standards are not provided with this test. As good laboratory practice, it is recommended to perform positive and negative controls periodically to verify test performance.

LIMITATIONS

This test is intended for the qualitative detection of SARS-CoV-2 virus antigen only. The exact concentration of SARS-CoV-2 viral antigen cannot be determined by this test.

Test results are for clinical reference only and should not be the sole basis for clinical diagnosis and treatment. Clinical management of patients should be considered in combination with their symptoms, physical signs, patient history, other laboratory tests, therapeutic responses and epidemiological information.

Proper sampling is crucial. Failure to follow the procedure can lead to incorrect test results. Improper collection, storage or even freezing and thawing of the sample can lead to inaccurate test results.

A false-negative test result may occur if the viral antigen level in a sample is below the detection limit of the test or if the sample was not collected or transported properly; therefore, a negative test result does not exclude the possibility of SARS-CoV-2 infection.

A positive result does not exclude co-infection with other pathogens.

Monoclonal antibodies may not detect SARS-CoV-2 viruses with slightly altered amino acid levels in the region of the target epitope, or may detect them with less sensitivity.

The amount of antigen in a sample may decrease with increasing disease duration. Samples collected after day 5 of illness are more likely to be negative compared to an RTPCR test.

The tests target the nucleocapsid proteins. Performance is not affected by mutations in the spike protein. Mutations in the nucleocapsid protein are not excluded in the future.

CHARACTERISTICS OF PERFORMANCE

The clinical performance of the *Green Spring® SARS-CoV-2 Antigen Rapid Test* was determined in prospective, randomised, single-blind studies. A total of 365 nasopharyngeal samples from symptomatic and asymptomatic patients were collected within 5 days of the onset of initial symptoms. The performance of the kit was compared with the results of a commercially available molecular test. The PCR comparisons use a nasopharyngeal swab.

Table 1: Clinical study nasopharyngeal (nose-throat)

Green Spring SARS-CoV-2 Antigen Rapid Test	PCR-Comparison		Total
	Positive	Negative	
Positive	150	0	150
Negative	5	210	215
Total	155	210	365
Sensitivity	96,77% (95%KI: 92,24-98,81%)		
Specificity	100,00% (95%KI: 97,76-100%)		
Accuracy	98,63% (95%KI: 96,89-100%)		

PPA(CI≤37): 96,77% (150/155), (95%KI: 92,24-98,81%)
NPA(CI≤37): 100,00% (210/210), (95%KI: 97,76-100%)

For the anterior nasal swab method, a total of 298 anterior nasal samples were collected from symptomatic and asymptomatic patients within 5 days of the onset of the first symptoms. The performance of

the kit was compared with the results of a commercially available molecular test. The PCR comparisons use a nasopharyngeal swab.

Table 2: clinical study anterior-nasal (nose-front)

Green Spring SARS-COV2 Antigen Rapid Test	PCR-Comparison		Total
	Positive	Negative	
Positive	154	0	154
Negative	6	138	144
Total	160	138	298
Sensitivity	96,25% (95%KI: 91,65-98,47%)		
Specificity	100,00% (95%KI: 96,62-100%)		
Accuracy	97,99% (95%KI: 96,97-100%)		

PPA(CI≤37): 96,25% (154/160), (95%KI: 91,65-98,47%)
 NPA(CI≤37): 100,00% (138/138), (95%KI: 96,62-100%)

For the saliva swab method, a total of 298 saliva samples from symptomatic and asymptomatic patients were collected within 5 days of the onset of the first symptoms. The performance of the kit was compared with the results of a commercially available molecular test. The PCR comparisons use a nasopharyngeal swab.

Table 3: clinical study saliva (lollipop)

Greenspring SARS-COV2 Antigen Rapid Test	PCR-Comparison		Total
	Positive	Negative	
Positive	147	0	147
Negative	13	138	151
Total	160	138	298
Sensitivity	91,88% (95%KI: 86,22-95,43%)		
Specificity	100,00% (95%KI: 96,62-100%)		
Accuracy	95,64% (95%KI: 93,32-97,96%)		

PPA(CI≤37): 91,88% (147/160), (95%KI: 86,22-95,43%)
 NPA(CI≤37): 100,00% (138/138), (95%KI: 96,62-100%)

		(Yes/No)
Influenza A	1.6 x 10 ⁵ TCID ₅₀ /mL	No
Influenza B	1.6 x 10 ⁵ TCID ₅₀ /mL	No
Human coronavirus HKU1	1.6 x 10 ⁵ TCID ₅₀ /mL	No
Human coronavirus OC43	1.6 x 10 ⁵ TCID ₅₀ /mL	No
Haemophilus influenzae	2.2x 10 ⁵ TCID ₅₀ /mL	No
MERS-coronavirus	2.1 x 10 ⁵ TCID ₅₀ /mL	No
SARS-coronavirus	3.2 x 10 ⁵ PFU/mL	Yes
Adenovirus C1	1.5 x 10 ⁵ TCID ₅₀ /mL	No
Adenovirus 71	1.5 x 10 ⁵ TCID ₅₀ /mL	No
Candida albicans	4.2 x 10 ⁵ CFU/mL	No
Respiratory syncytial virus	5.1 x 10 ⁵ TCID ₅₀ /mL	No
Enterovirus	5.4 x 10 ⁵ TCID ₅₀ /mL	No
Malaria	2.2 x 10 ⁶ CFU/mL	No
Dengue	1.2 x 10 ⁵ TCID ₅₀ /mL	No
Human coronavirus NL63	1.7x 10 ⁵ TCID ₅₀ /mL	No
Human coronavirus 229E	2.2 x 10 ⁵ TCID ₅₀ /mL	No
Streptococcus pneumoniae	1.1 x 10 ⁶ CFU/mL	No
Pneumocystis jirovecii (PJP)	1.0 x 10 ⁵ TCID ₅₀ /mL	No
Legionella pneumophila	1.4 x 10 ⁶ CFU/mL	No
Chlamydia pneumoniae	1.1 x 10 ⁶ IFU/mL	No
Human Metapneumovirus (hMPV)	1.1 x 10 ⁵ TCID ₅₀ /mL	No
Parainfluenza virus 1	1.0 x 10 ⁵ TCID ₅₀ /mL	No
Parainfluenza virus 2	1.0 x 10 ⁵ TCID ₅₀ /mL	No
Parainfluenza virus 3	3.5 x 10 ⁵ TCID ₅₀ /mL	No
Parainfluenza virus 4	1.4 x 10 ⁵ TCID ₅₀ /mL	No
Rhinovirus	1.3 x 10 ⁵ PFU/mL	No
Mycoplasma pneumoniae	1.8 x 10 ⁶ CFU/mL	No
Bordetella pertussis	1.5 x 10 ⁶ CFU/mL	No
Mycobacterium tuberculosis	1.0 x 10 ⁶ CFU/mL	No
Concentrated human nasal contents (representative of normal respiratory microbial flora)	100%	No
Streptococcus pyogenes	1.0 x 10 ⁶ CFU/mL	No



CROSS-REACTIVITY

Potential cross-reactant	Concentration	Cross-reactivity
--------------------------	---------------	------------------

INTERFERENCE

SARS-CoV-2 antigen nasal swab samples were mixed with one of the following substances to specific concentrations and tested in multiple replicates. No false-positives or false-negatives were found:

Substance	Concentration	Substance	Concentration
Whole Blood	5%	Naso GEL(Nei Med)	6%v/v
Fluticasone Propionate	4%v/v	Mucin	0.54%
CVS Nasal Drops(Phenylephrine)	17%v/v	Ricola(Menthol)	1.6mg/mL
Tamiflu (Oseltamivir Phosphate)	6mg/ml	Afrin (Oxymetazoline)	14%v/v
Sucrets (Dyclonin/Menthol)	1.4 mg/mL	CVC Nasal Spray(Cromolyn)	16%v/v
Chloraseptic (Menthol/Benzocaine)	1.8 mg/mL	Nasal Gel (Oxymetazoline)	9%v/v
Homeopathic(Alkalol)	1:10dilution	Mupirocin	12 mg/mL
Ore Throat Phenol Spray	16%v/v	Fisherman's Friend	1.3mg/mL
Tobramycin	5 µg/mL	Zicam	4%v/v

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In-vitro diagnostic use	Read the instruction before using	CE marking
Batch number	Expiry date	Date of manufacture
Do not re-use it	Store at 2C -30C	Keep away from sunlight
Keep dry	Manufacturer	EU authorised representative

LIMIT OF DETECTION (ANALYTICAL SENSITIVITY)

The limit of detection (LOD) for the *Green Spring® SARS-CoV-2 Antigen Rapid Test* is 4×10^2 TCID₅₀/mL. The LOD for *Green Spring® SARS-CoV-2 Antigen Rapid Test Kit* was determined using limiting dilution of a gamma irradiation inactivated virus sample. The sample was provided at a concentration of 1.3×10^6 TCID₅₀/mL.

HIGH-DOSE HOOK EFFECT

The LOD study tested the highest concentration of the sample (TCID₅₀ of 1.3×10^6 TCID₅₀/mL). No hook effect was observed.

FURTHER PRODUCT INFORMATION

Manufacturer: Shenzhen Lvshiyuan Biotechnology Co., Ltd

101,201,301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New District, Shenzhen, 518120 China

EU representative: Obelis s.a.

Bd General Wahis 53, 1030 Brussels Belgium

Importer: Better AG

General-Guisan-Str. 8
6300 Zug, Switzerland
Tel: + 353 1 513 7511
Email: info@OdemShop.com
Shop: www.OdemShop.com

Authorization

It is hereby certified and declared that the company:

Our Shenzhen Lvshiyuan Biotechnology Co.,Ltd. is a legally established and

-manufacturing(COVID- 19(2019-nCoV) Coronavirus IgG/IgM Rapid Test Kit、

SARS-CoV-2 Antigen Rapid Test Kit 、ARS-CoV-2 Neutralizing Antibody Rapid Test Kit)

based enterprise Its main place of business is located in 101,201,301,D Building,

“Better AG” located in General-Guisan-Str. 8, 6300 Zug, Switzerland Is authorized to import, sell, distribute the "Green spring" branded goods in Europe, Asia and Africa.

We hereby confirm the authenticity of the antigene tests sold by this distributor.

The authorization period is from May 24,2021 to May 23.2025

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Shenzhen Lvshiyuan Biotechnology Co., Ltd

Add: D Building, National Biological Industrial Park of Marinelife, No.2 Binhai Road, Dapeng, Shenzhen, 518120, China.

The Statement on detection of mutant viruses

WHO held an emergency meeting on 26 November 2021 to discuss the recently discovered mutant strain of the novel coronavirus B.1.1.529. After the meeting, WHO issued a statement, has designated B.1.1.529 as a "Variant of Concern" , named Omicron.

The B.1.1.529 variant was first reported to WHO from South Africa on 24 November 2021 and the first sample infected with the mutant strain was collected on 9 November, the WHO said in a statement .This variant has a large number of mutations, some of which are concerning

Preliminary studies suggest that this variant causes an increased risk of reinfection in humans compared to other "concerns" variants. The number of cases of this variant appears to be increasing in almost all provinces in South Africa.

Since the outbreak of the SARS-CoV-2 virus, it has been reliably reported that there have been at least hundreds of mutations in the gene sequence, all of which have resulted in the virus being more infectious and more pathogenic. The most famous strains include Alpha, which first appeared in Kent, England, Beta in South Africa, Gamma in Brazil, Delta in India and Mu which was first discovered in Colombia and officially named by the WHO on August 30. And Omicron, recently discovered in South Africa. Among them, Omicron variant virus has recently invaded many countries and regions around the world with its strong infectivity and pathogenicity.

The new coronavirus (SARS-CoV-2 or 2019-nCoV) is a non-segmented forward RNA virus. This is the cause of the new type of coronavirus pneumonia (COVID-19), which is highly contagious in humans. The SARS-CoV-2 virus has several structural proteins, including spikes (S), envelope (E), membrane (M) and nucleocapsid (N).

The SARS-CoV-2 virus has the characteristics of strong nucleocapsid (N) protein stability. The mutant virus strains that have been found worldwide are derived from the SARS-CoV-2 20B/GR evolutionary strain (lineage B.1.1.7), including many mutation, the mutation location is the spike (S) protein of the new coronavirus, which is the location where the SARS-CoV-2 virus uses to bind to the cell's ACE2 receptor.

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Shenzhen Lvshiyuan Biotechnology Co., Ltd

Add: D Building, National Biological Industrial Park of Marinelife, No.2 Binhai Road, Dapeng, Shenzhen, 518120, China.

The SARS-CoV-2 Antigen Rapid Test Kit produced by Shenzhen Lvshiyuan Biotechnology Co., Ltd. is used for in vitro qualitative detection of SARS-CoV-2 virus nucleocapsid (N) protein in human nasopharyngeal , oropharyngeal , anterior -nasal or saliva samples.

It can be seen that the mutation sites of mutated virus strains including Omicron strain have no effect on the detection rate of the kits produced by our company. The kit is suitable for assay of the SARS-CoV-2 variant virus called 'Omicron '.

Shenzhen Lvshiyuan Biotechnology .Co., Ltd.

Chief Director :Mr Jiang

Date: 28th November .2021

Tel: +86-755-28438788

Fax: +86-755-28938800

ISO 13485 认证证书 ISO 13485 certification

SGS

Certificate CN20/42084

The management system of

Shenzhen Lvshiyuan Biotechnology Co., Ltd.

101, 201, 301, Building D, No. 2, Industrial Avenue, Buxin Village,
Buxin Community, Dapeng Subdistrict Office, Dapeng New District,
Shenzhen, Guangdong, 518120, P.R. China

has been assessed and certified as meeting the requirements of

ISO 13485:2016 EN ISO 13485:2016

For the following activities

**Design, Manufacture and Distribution of Dry Fluorescent
Immunoassay Instrument and In Vitro Diagnostic Test Kits
(ELISA, Colloidal Gold) for SARS-CoV-2, Influenza A and Influenza B.**

This certificate is valid from 2 March 2021 until 13 June 2023
and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 26 May 2023
Issue 2. Certified since 14 June 2020

Authorised by

SGS United Kingdom Ltd
Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK
t +44 (0)151 350-6666 f +44 (0)151 350-6600 www.sgs.com

HC SGS 13485 2016 0118

Page 1 of 1

0005

This document is issued by the Company subject to its General Conditions of Certification Services accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <http://www.sgs.com/certified-clients-and-products/certified-clients-directory>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

ISO9001 认证证书 ISO9001 certification

CERTIFICATE OF REGISTRATION

The Quality Management Systems of

Shenzhen Lvshiyuan Biotechnology Co., Ltd

Unified Social Credit Code:914403007576264357

Registration address:101, 201, 301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New District, Shenzhen

Production address:D Building, National Biological Industrial Park Of Marinelife, Binhai No.2 Road, Dapeng, Shenzhen

has been assessed by GIC and complying with

GB/T19001-2016/ISO9001:2015

For the following activities

**Research and development, production and service of
food safety testing kits, animal disease diagnostic kits
and test cards**

Date of Issue: 13 February 2019

Date of Expiry: 12 February 2022

Date of Initial Certification: 13 February 2019

Certificate No.: J19Q2GZ8012523R0M

Scan for certificate status

The granting of this certificate does not mean that the certificate holder can avoid any legal obligation. If the products or activities covered in the scope of certification require administrative license, the certificate shall be only valid within the scope of administrative licensing. The registered organization shall be subject to regular annual supervision by GIC, and the continual validity of the certificate is base upon conformity of audit. Please scan two-dimension code at left to find the certificate information. This certificate can be queried at Certification and Accreditation Administration of the People's Republic of China official website (www.cnca.gov.cn) & GIC website (www.gicg.com.cn)

GIC WeChat public number

Signature:

Guardian Independent Certification Ltd

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Germany BfArM Approval of antigen tests for professional testing

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12.02.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with $CT < 25$, 23 samples with CT between 25 and 30, and 9 samples with $CT > 30$. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 μ L of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that many other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

Contact

Email: sarscov2ivd@pei

Last updated: 12.02.2021

Overview of SARS-CoV-2 Antigen Rapid Tests Assessed as Reflecting the Current State of the Art

Name of Test	Manufacturer (Distributor)
Panbio™ COVID-19 Ag Rapid Test Device (NASOPHARYNGEAL)	Abbott Rapid Diagnostics Jena GmbH
RIDA®QUICK SARS-CoV-2 Antigen	R-Biopharm AG
SARS-CoV-2 Rapid Antigen Test	SD BIOSENSOR (Roche Diagnostics GmbH)
NADAL® COVID-19 Ag Schnelltest	nal von minden gmbh
STANDARD™ F COVID-19 Ag FIA	SD BIOSENSOR
STANDARD™ Q COVID-19 Ag Test	SD BIOSENSOR
BIOSYNEX COVID-19 Ag BSS	BIOSYNEX SWISS SA
MEDsan® SARS-CoV-2 Antigen Rapid Test	MEDsan GmbH
TestNOW® - COVID-19 Antigen	Affimedix
NowCheck® COVID-19 Ag Test	BIONOTE
Coronavirus Ag Rapid Test Cassette (Swab)	Zhejiang Orient Gene Biotech Co.,Ltd
Sofia SARS Antigen FIA	Quidel Corporation
COVID-19 Ag Test Kit	Guangdong Wesail Biotech Co., Ltd.
CLINITEST® Rapid COVID-19 Antigen Test	Siemens Healthineers
ESPLINE® SARS-CoV-2	Fujirebio Inc. (Mast Diagnostica GmbH)
BD Veritor™ System for Rapid Detection of SARS-CoV-2	Becton Dickinson
GenBody COVID-19 Ag	IVC Pragen Healthcare
LumiraDx SARS-CoV-2 Ag Test	LumiraDX
Exdia COVID-19-Ag-Test	Precision Biosensor Inc. (Axon Lab AG)
SARS-CoV-2 Ag Rapid Test (FIA)	Wantai (Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.)
SARS-CoV-2 Antigen Schnelltest	Xiamen Boson Biotech Co., Ltd (Medicovid-AG; technomed GmbH; Löwe Medizintechnik)
COVID-19 Antigen Schnelltest (Colloidal Gold)	Joinstar Biomedical Technology Co., Ltd (CIV care impuls Vertrieb)
mö-screen Corona Antigen Test	Mölab GmbH
Rapid SARS-CoV-2 Antigen Test Card	MP Biomedicals Germany GmbH
Lyher Novel Coronavirus (COVID-19) Antigen Test Kit (Colloidal Gold)	Hangzhou Laihe Biotech Co., Ltd. (Lissner Qi GmbH)
AMP Rapid Test SARS-CoV-2 Ag	Ameda Labordiagnostik GmbH
Clungene COVID-19 Antigen Rapid Test	Hangzhou Clongene Biotech Co., Ltd.
Gensure™ COVID-19 Antigen Rapid Test Kit	GenSure Biotech Inc.
SARS-CoV-2 Antigen Rapid Test Kit	Beijing Lepu Medical Technology Co., Ltd
Hightop SARS-CoV-2 (Covid-19) Antigen Rapid Test	Qingdao Hightop Biotech Co., Ltd.
Rapid Covid-19 Antigen Test (Colloidal Gold)	Anbio (Xiamen) Biotechnology Co., Ltd

Name of Test	Manufacturer (Distributor)
Safecare COVID-19 Ag Rapid Test Kit (Swab)	Safecare Biotech Hangzhou Co., Ltd.
QuickProfile Covid-19 Antigen Test Card	LumiQuick Diagnostics, Inc.
Covid 19 Antigen Schnelltest	BioRepair GmbH
Green Spring SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	Shenzhen Lvshiyuan Biotechnology Co., Ltd.
CAT Antigen Covid Rapid Test	Oncosem Onkolojik Sistemler San. Ve Tic. A.S.
ScheBo SARS-CoV-2 Quick Antigen	ScheBo Biotech AG
Nova Test SARS-CoV-2 Antigen Rapid Test Kit	Atlas Link Technology Co., Ltd.
Toda Coronadiag Ag	Toda Pharma
Humasis COVID-19 Ag Test	Humasis Co., Ltd.
Beijing Hotgen Biotech Co., Ltd.	Neuartiges Coronavirus (2019-nCoV)-Antigentest (Kolloidales Gold); Novel Coronavirus 2019-nCoV Antigen Test (Colloidal gold)
Xiamen AmonMed Biotechnology Co.,Ltd.	COVID-19 Antigen Rapid Test Kit (Colloidal Gold)
Canea COVID-19 Antigen Schnelltest	Core Technology Co., Ltd.
fluorecare COVID-19 SARS-CoV-2 Spike Protein Test Kit (Colloidal Gold Chromatographic Immunoassay)	Shenzhen Microprofit Biotech Co., Ltd
Tetsealabs® Rapid Test Kit COVID-19 Antigen Test Cassette	Hangzhou Testsea Biotechnology Co., Ltd
Lysun COVID-19 Antigen Rapid Test Device (Colloidal Gold)	Hangzhou Lysun Biotechnology Co., Ltd.

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Obtained screenshots of foreign standard certification

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欧盟通用和互认清单 On the EU common list & Mutual Recognition list

https://covid-19-diagnostics.jrc.ec.europa.eu/devices/detail/2109

COVID-19 In Vitro Diagnostic Devices and Test Methods Database

Home > COVID-19 In Vitro Diagnostic Medical Devices > COVID-19 In Vitro Diagnostic Medical Device - detail

COVID-19 In Vitro Diagnostic Medical Device - detail

< Previous Next >

Green Spring SARS-CoV-2 Antigen-Rapid test-Set

Manufactured by Shenzhen Lvshiyuan Biotechnology Co., Ltd., China - <https://www.lsybt.com/>

Device identification number	2109
CE Marking	Yes
HSC common list	Yes
HSC mutual recognition	Yes
Format	Manual, Near POC / POC
Physical Support	Lateral flow
Target	Antigen
Specimen	Anterior nasal swab, Nasal swab, Nasopharyngeal swab, Oropharyngeal swab
Pathogens detected	Coronaviruses (HCoV), SARS-CoV
Lineages detected	B.1.1.7 (United Kingdom), B.1.351 (South Africa), P.1 (Japan/Brazil)
Commercial Status	Commercialised
Last Update	2021-07-07 05:24:52 CET
Comments	Der Green Spring® SARS-CoV-2-Antigen-Schnelltest dient dem schnellen qualitativen Nachweis des Nukleocapsid-Protein-Antigens von SARS-CoV-2 in menschlichen Nasen-, Nasen-Rachen oder Rachenabstrichproben. Die Ergebnisse dienen dem Nachweis von SARS-CoV-2-Antigenen.

比利时白名单 Whitelist of Belgium

Manufacturer	Device Name	CE Marking	HSC Common List	HSC Mutual Recognition	Format	Physical Support	Target	Specimen	Pathogens detected	Lineages detected	Commercial Status	Last Update	Comments	Whitelist Status
LiClear Biotech (Hangzhou)	SARS-CoV-2 Nucleocapsid (N) Antigen Rapid Test Cassette	95.1	99.3					NP swab	-					No
LumiQuick Diagnostics	QuickProfile COVID-19 Antigen Test Strip	94.0	99.0					NP swab						Yes
MEDaan	SARS-CoV-2 Antigen Rapid Test	92.9	99.8					Nasal swab/OP swab						Yes
MP Biomedicals	Rapid SARS-CoV-2 Antigen Test Card	96.4	99.0					NP swab/OP swab						Yes
Multi-G	COVID19CHECK-SEN	92.0	99.2					NP swab/OP swab						No
	COVID19CHECK-NAS (COVID-19 antigen rapid test cassette)	96.1	98.2					Nasal swab						No
	COVID19CHECK-SAL (Sars-Cov-2 Antigen Rapid Test)	99.8 (Cis 33)	99.1					Saliva						No
MyLab Discovery Solutions Pvt	PathoCatch COVID-19 Antigen Lateral Flow Test Device	92.0	100.0					Nasal swab						No
hel von minden	bioDIAL COVID-19 Ag Rapid Test	97.8	99.9					NP swab/OP swab						Yes
Nanjing Vazyme Medical Technology	Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) Antigen Detection Kit (Colloidal Gold-Based)	97.8	99.2					Nasal swab/OP swab						18.8.9
NanoCheck	COVID-19 Antigen Test Kit (Colloidal Gold) (Cassette)	94.7	100.0					NP swab						14.10
Nantong Diagnosis Biotechnology	COVID-19 Antigen Test Kit (Colloidal Gold)	96.1	99.3					Saliva						16.1.7
New Gene (Hangzhou) Biotechnology	COVID-19 Antigen Test Kit (Colloidal Gold)	98.6/98.5	100.0					NP swab/OP swab						No
ONCOSEM Onkoloogik Szisztémés és Tit. Á.S.	COVID-19 Antigen Test Kit	97.3/95.7/95.1	99/99/99.1					Nasal swab/OP swab/Sputum						15.01
Onthe-Clinical Diagnostics	2019-nCoV Antigen Rapid Test Kit	97.6	99.3					NP swab						-
PCI	WITROS SARS-CoV-2 Antigen test	97.8	99.2					NP swab						Lab test
PCI	PCI COVID19 Ag Gold	80.8 (CIS 30)/91.7	99.5/100					Saliva/NP swab						17.10
Prestige Diagnostics	2019-nCoV Antigen Device	96.9	98.1					NP swab						13.10
PRIMA Lab	PRIMACOV COVID-19 Antigen Rapid Test	95.4/92.8	99.2/100					NP swab/Nasal swab						17.8.7
Qimedao Hiitop Biotech	SARS-CoV-2 Antigen Rapid Test (Immunochromatography)	92.7/95.0/95.0	99.8/99.8/99.8					Nasal swab/OP swab/OP swab						13.11
Quidel	Sofia SARS-CoV-2 Antigen RIA	96.7	100.0					NP swab/Nasal swab						10.8.7
	Sofia 2 Plus - SARS Antigen P1	95.2	100.0					NP swab/Nasal swab						14.10
Roche Diagnostics	Elecsa SARS-CoV-2 Antigen	94.5	99.9					NP swab/OP swab						14.10
SD Biosensor (distributed by Roche)	SARS-CoV-2 Rapid Antigen Test	96.5	99.7					NP swab						16.10
Shenzhen Huan Bios Technology	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	94.6	99.1					Nasal swab						No
Shenzhen Landwind Biotechnology	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	95.7	99.2					Nasal swab						No
Shenzhen Lvshiyuan Biotechnology	Green Spring SARS-CoV-2 Antigen-Rapid Test Kit (Colloidal Gold)	96.6/99	100.0					Nasal swab/OP swab						21.09
Shenzhen Mikropoint Biotech	Fluoresce SARS-CoV-2 Test Kit	98.0/100	100.0					NP swab/Nasal swab						No
	Fluoresce SARS-CoV-2 Spike Protein Test Kit	96.9	100.0					NP swab						12.10
Shenzhen Ultra Diagnostics Biotech	SARS-CoV-2 Antigen Test Kit (SCC001 + SCC002)	92.100/98	100.0					NP swab/OP swab/Saliva						No
	SARS-CoV-2 Antigen Test Kit (SCC003 + SCC004)	92/96.7/97.2	100/99/99					NP swab/OP swab/Sputum						No
Shenzhen YHLO Biotech Co	GLINE-2018-nCoV Ag	95.8/92.8	99.3/99.3					NP swab/Nasal swab						13.17
SinScreen	COVID-19 Coronavirus Rapid Antigen Test Cassette	93.2	100.0					NP swab/OP swab						17.7.5
Tedapharma	TODA Coronadiag Ag	96.6	100.0					NP swab/OP swab						14.10
Ulf med Products	COVID-19 Antigen Test (Nasopharyngeal Swab)	96.4	99.2					NP swab						13.10
	COVID-19 and Influenza A/B Antigen Combo Rapid Test (Nasopharyngeal Swab)	96.4	99.2					NP swab						13.10
Van Oosteen Medical	Coronavirus Ag Rapid Test Cassette (Swab)	96.7	98.2					NP swab						16.11
VivaDiag Pro	VivaDiag Pro SARS-CoV-2 Ag Rapid Test	96.1/96.1/97.0	100.0					NP swab/OP swab/Nasal swab						21.03
VivaCheck Biotech	VivaDiag SARS-CoV-2 Ag Rapid Test	95.0	100.0					NP swab/OP swab/Nasal swab						21.02
	VivaDiag SARS-CoV-2 Ag Saliva Rapid Test	96.3	100.0					Saliva						18.10

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
英国注册证书 Registration Certificate of UK

Medicines & Healthcare products
Regulatory Agency

Medicines & Healthcare products
Regulatory Agency

10 South Colonnade
Canary Wharf
London
E14 4PU
United Kingdom

gov.uk/mhra

PureUKCA Ltd
59
St. Martin's Lane
Middlesex
London
WC2N 4JS
England, United Kingdom

09 September 2021

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We are pleased to confirm that the application to register or update an existing registration for the following manufacturer, which you submitted on **09 September 2021** has been reviewed:

Application reference: **2021090901215399**

Manufacturer organisation: **Shenzhen Lvshiyuan Biotechnology Co.,Ltd.**

Address:

**D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng
New District, Shenzhen, 518120 China
Shenzhen
518120
China**

Manufacturer registration status: **Registered**

Device(s):

GMDN term	Status	Comment
SARS-CoV-2 antigen IVD, kit, immunochromatographic test (ICT), rapid	Registered	

Please note this letter **does not** represent any form of accreditation, certification or approval by the UK Competent Authority.

If you stop placing devices on the market or if you are not complying with the Regulations, you should inform us so that we can amend our records. You should be aware that it is an offence to place on the market UKCA or CE marked devices that do not comply with the regulations.

Please inform us of the following chargeable changes:

- 1. company/organisation information e.g. name and address**
- 2. additional devices (GMDN code or term)**

Please also use the Devices Online Registration Database (DORS) to tell us of the following changes e.g. removal/discontinuation of a device (GMDN) or product from your registration record, change of contact person, telephone number and/or email address, for which payment of our statutory fee does not apply.

Please note that the name and address of manufacturer, UK Responsible Person or Authorised Representative (Northern Ireland only) and devices that have been registered will be published on our [Public Access Registration Database \(PAR\)](#).

The account number for your company/organisation is **0000018481**.

Yours sincerely,

Ngozi Onyeukwu
Device registrations service
Devices division
MHRA

荷兰注册证 Registration Certificate in Netherlands

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
Ref No.: BS 0171-2020

Annex A - List of Devices

(Recital 29 of the Directive 98/79/EC on In Vitro Diagnostic Medical Devices)

#	Catalogue reference number	Commercial Name	Generic Device Term	Short description and intended use	GMDN/EDMS Code	Class
1	GF102B1	SARS-CoV-2 Antigen Rapid Test Kit	SARS-CoV-2 (COVID-19) Antigen Rapid Test Kit (Colloidal Gold)	SARS-CoV-2 (Covid-19) Antigen Rapid Test Kit is used to quantitatively detect antigen of the novel coronavirus in nasal swab.	15.04.80.19	Others

CE 认证 CE certification

CERTIFICATE OF IVD NOTIFICATION

Ref. No.: BS 0171-2020

BELGIUM

Date: 19/11/2020

Order No.: OG 0117-2020

THIS IS TO CERTIFY THAT, ACCORDING TO THE COUNCIL DIRECTIVE 98/79/EC, OBELIS S.A. (O.E.A.R.C.) PERFORMED ALL NOTIFICATION DUTIES AND RESPONSIBILITIES AS THE EUROPEAN AUTHORIZED REPRESENTATIVE (EC REP) OF:

NAME: SHENZHEN LVSHIYUAN BIOTECHNOLOGY CO., LTD.

ADDRESS: 101, 201, 301, D BUILDING, NO. 2 INDUSTRIAL AVENUE, BUXIN VILLAGE, BUXIN COMMUNITY, DAPENG SUBDISTRICT OFFICE, DAPENG NEW DISTRICT, SHENZHEN, 518120, CHINA

AS STIPULATED AND DEMANDED BY THE AFOREMENTIONED DIRECTIVE.

The Manufacturer declares that the IVD devices comply with the Directive including all essential requirements.

The Manufacturer has provided Obelis s.a. (O.E.A.R.C.) with all the appropriate declarations according to the 98/79/EC Directive – article 10 requirements including the EC Declaration of Conformity confirming that his In-Vitro Diagnostics medical devices, as stipulated here above, are fulfilling the applicable requirements of the European Council Directive 98/79/EC

The notification of the following In-Vitro Diagnostic medical devices has been completed by Obelis s.a. (O.E.A.R.C.) on the 18/11/2020 in compliance with the European Council Directive 98/79/EC - article 10 requirements.

IN-VITRO DIAGNOSTIC MEDICAL DEVICES: PLEASE SEE ANNEX A - LIST OF DEVICES (1 PAGE, 1 DEVICE)

As of the 19/11/2020, and as long as the manufacturer will continue complying with the hereabove mentioned requirements* he therefore:

- Is required to affix the CE marking on these devices;
- Place these devices in the Territory of Belgium and/or the other EEA Member States (excluding territories not in alignment with Decision 2010/227/EU).

Obelis s.a. - O.E.A.R.C.
Registered Address:
Rue de la Woluwe 62/63
1200 Brussels
Belgium

Mr. G. Elkayam CEO

Obelis sa

Obelis European Authorized Representative Center is a member of the European Association of Authorized Representatives (E.A.A.R.), ISO 9001 : 2015 and ISO 13485 : 2016 certified in accordance to the profession of a European Authorized Representative.

* This is not a CE mark and is only provided as a template for informational purposes.

** This Certificate will be automatically void if the notification is rejected by the EU Authorities or upon termination of the EAR agreement.

Registered Address : Bcd. Général Wahnis 53-1090 Brussels | Registered Office Address: Bd Brand Whitlock 30, B-1200 Brussels - Belgium

Declaration of Conformity

According to annex III of the Council Directive 98/79/EC on in vitro diagnostic medical device
We,

Company Name: Shenzhen Lvshiyuan Biotechnology Co., Ltd.

*Address: 101, 201, 301, D Building, No. 2 Industrial Avenue, Buxin Village, Buxin Community,
Dapeng Subdistrict Office, Dapeng New District, Shenzhen 518120 China*

Declare under our sole responsibility that the following in vitro diagnostic medical devices
other than those covered by annex II and devices for performance evaluation

List of Products:

1. SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Meet the provisions of the Council Directive 98/79/EC concerning medical devices which
apply to them.

Undersigned declares to fulfill the obligations imposed by Annex III section 2 to 5:

- availability of the technical documentation set in Annex III (section 3), allowing the assessment of conformity of the product with the requirements of the Directive.
- the manufacturer shall take necessary measures to ensure that the manufacturing process follows the principles of quality assurance as appropriate for the products manufactured (Annex III section 4).
- the manufacturer shall institute and keep up to date a systematic procedure to review experience gained from devices in the post-production phase and to implement appropriate means to apply any necessary corrective actions (Annex III section 5).

Conformity assessment was performed according to Article 9 (7) and Annex III, section 3.

Our current Quality System is formatted to international standards:

- ISO 9001: 2015

Corporate Contact Information

COMPANY NAME: Shenzhen Lvshiyuan Biotechnology Co., Ltd.

COMPANY ADDRESS: 101, 201, 301, D Building, No. 2 Industrial Avenue, Buxin Village,

RESPONSIBLE PERSON'S name: Jiang Yongqing

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依据联合国 GHS 制度第八修订版编/According to UN GHS (the 8th revised edition)

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ISSUED BY: ZHONGKE SERVICES OF TESTING (SHENZHEN) CO., LTD.


2022 年 01 月 印制 Jan. 2022 Print

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BETTER AG

General-Guisan-Str. 8
6300 Zug, Switzerland

Contact:

EN: + 353 1 513 75 11
CH: + 41 (0) 71 58 80 248
www.OdemShop.com
E-Mail: info@OdemShop.com